

## STANDARD WHEELCHAIRS

# ADVANTAGE® LX WHEELCHAIR

### LX SERIES FIXED BACK WHEELCHAIRS

The Advantage® Wheelchair provides comfort and durability at an affordable price. Designed as a low maintenance wheelchair, it offers an attractive silver vein durable hammertone steel finish.

#### Features


- Available in 16", 18" and 20" seat widths
- Black padded nylon and vinyl upholstery available with standard chart pocket
- Carbon steel frame provides strength and durability
- Industry compatible tool-free front rigging with heel loops standard on swingaway footrests
- Hammertone paint provides a long lasting and durable finish
- Comfortable padded armrests
- Maximum Weight Capacity 300 lb – EVENLY DISTRIBUTED
- Meets Medicare code K0001

#### Warranty:

One (1) year limited warranty



*fixed-back model with detachable desk arms (shown) and elevating legrest (shown) – other configurations available*



## Advantage® LX Series Fixed-Back Wheelchair

Item	HCPCS+ code	Overall Chair Width		Overall Chair Depth		Overall Chair Height		Folded Chair Width		Seat Width		Seat Depth		Seat Height		Back Height		Wheel Diameter		Caster Diameter		Product <sup>®</sup> Weight w/o front rigging		Weight <sup>†</sup> Capacity EVENLY DISTRIBUTED		Fixed Full	Detachable Desk	Detachable Full	Swingaway Footrest	Elevating Legrest	Nylon	Vinyl	Rear Anti-Tippers
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kgf	lb	kgf	arms	front rigging	upholstery					
3H020100	+E0951	26.5	67.3	32	81.3	36	91.4	11.5	29.2	18	45.7	16	40.6	19.5	49.5	16.5	41.9	24	61.	8	20.3	40	18.1	300	136.1	■	—	—	■	—	■	—	■
3H020110	+E0990 OR K0195																									■	—	—	■	—	■	—	■
3H020120	+E0951																									—	■	—	■	—	■	—	■
3H020130	+E0990 OR K0195																									—	■	—	■	—	■	—	■
3H020140	+E0951																									—	■	—	■	—	■	—	■
3H020220	+E0951	26.5	67.3	32	81.3	36	91.4	11.5	29.2	16	40.6	16	40.6	19.5	49.5	16.5	41.9	24	61.	8	20.3	40	18.1	300	136.1	—	■	—	■	—	■	—	■
3H020230	+E0990 OR K0195																									—	■	—	■	—	■	—	■
3H020320	+E2201+E0951	30.5	77.5	32	81.3	36	91.4	11.5	29.2	20	50.8	16	40.6	19.5	49.5	16.5	41.9	24	61.	8	20.3	40.2	18.2	300	136.1	—	■	—	■	—	■	—	■
3H020330	+E2201+E0990 OR K0195																									—	■	—	■	—	■	—	■
3H020340	+E2201+E0951																									—	■	—	■	—	■	—	■
3H020350	+E2201+E0990 OR K0195																									—	■	—	■	—	■	—	■
3H021100	+E0951	26.5	67.3	32	81.3	36	91.4	11.5	29.2	18	45.7	16	40.6	19.5	49.5	16.5	41.9	24	61.	8	20.3	40	18.1	300	136.1	■	—	—	■	—	■	—	■
3H021110	+E0990 OR K0195																									■	—	—	■	—	■	—	■
3H021120	+E0951																									—	■	—	■	—	■	—	■
3H021130	+E0990 OR K0195																									—	■	—	■	—	■	—	■
3H021220	+E0951	26.5	67.3	32	81.3	36	91.4	11.5	29.2	16	40.6	16	40.6	19.5	49.5	16.5	41.9	24	61.	8	20.3	40	18.1	300	136.1	—	■	—	■	—	■	—	■
3H021230	+E0990 OR K0195																									—	■	—	■	—	■	—	■
3H021320	+E2201+E0951	30.5	77.5	32	81.3	36	91.4	11.5	29.2	20	50.8	16	40.6	19.5	49.5	16.5	41.9	24	61.	8	20.3	40.2	18.2	300	136.1	—	■	—	■	—	■	—	■
3H021330	+E2201+E0990 OR K0195																									—	■	—	■	—	■	—	■

■ Standard ■ Optional ■ Available ■ Accessory — Not Apply / Available

\* designed to meet the requirements of † standard rounding method used † maximum weight capacity - evenly distributed - refer to user manual  
 § approximate measurement; product weight — as seat width and rigging change, so does product weight ★ options available at additional cost - accessories sold separately

### Parts & Accessories: ★ sold separately

	HCPCS+	UOM	Item No.
Footrest Heel Loop	E0951	pair	90763064
Wheel Lock Fixed Arm (with rubber tips)	—	pair	90763070
Wheel Lock Lever Extension	E0961	pair	90763062

	HCPCS+	UOM	Item No.
Positioning Belt (auto-style, push button) <small>(adjustable length: 27 - 47 inches / 68.58 - 119.38 cm)</small>	E0978	each	300-505
Combination O2 & IV Holder	E2208	each	B10320F
O2 Cylinder Holder	E2208	each	B10325
O2 Holder Bag (nylon)	E2208	each	EJ02BAG07

	HCPCS+	UOM	Item No.
Anti-tippers, detachable-arm sideframe	E0971	pair	90763565
Anti-tippers, fixed-arm sideframe	E0971	pair	90763566
Anti-Fold Bar (Chromed) Lower Frame	—	each	92523850
Wheelchair Positioning Tray	—	each	EJTRAY18



GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2019, GF Health Products, Inc. All Rights Reserved. Everest & Jennings is a registered trademark of GF Health Products, Inc.



GF1700099RevC24