

## Lab Chair

### MODEL 2166W

- Contemporary styling
- Cast aluminum base with satin finish
- 2" thick padded seat and backrest
- Pneumatic cylinder height adjustment
- Fully adjustable tilting seats
- Adjustable backrest with high impact back shroud
- Soft roll, dual wheel, hooded casters
- Adjustable foot ring
- Premium, stain-resistant, woven, knit-backed, vinyl upholstery
- 250 lbs. load capacity under normal use



California residents please be advised of the following, pursuant to Proposition 65:



**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

#### SPECIFICATIONS

Seat Diameter	18"
Base Diameter	23"
Height Range	21"-28.5"

#### Options



23 Optional Glides B

#### Upholstery Colors

