

Techno-Aide's one-of-a-kind **EZ Fit Sizing Tool** quickly and accurately fits people of all shapes and sizes into the perfect Techno-Aide apron size without the need for measuring tapes. However, use the provided measuring charts to verify your fit if a measuring tape is the only available option.


**Techno-Aide standard aprons (front, full wrap, vests, & skirts) are fit using two measurements:**

- Width measured at *the widest point of the body* applicable to the apron style
- Length measured applicable to the apron style




Visit [Techno-Aide.com/ez-fit-sizing-tool](https://Techno-Aide.com/ez-fit-sizing-tool) for step-by-step measuring instructions and to order your sizing tool today!

## STANDARD APRONS (FRONT, FULL WRAP, VESTS, & SKIRTS)

### WIDTHS (FRONT, FULL WRAP, VESTS, & SKIRTS)



### LENGTHS (FRONT, FULL WRAP, VESTS, & SKIRTS)





**SPECIALTY APPAREL**

**THYROID GUARDS**

MEASUREMENT	STANDARD	MAMMOGRAPHY	PEDIATRIC
LENGTH <i>end to end</i>	24"	24"	16"
BIB HEIGHT	6"	4"	4"

**PEDIATRIC FRONT APRONS**

MEASUREMENT	STANDARD
WIDTH <i>at widest point</i>	18"
LENGTH <i>from shoulders</i>	25.5"

**BLANKET GUARDS**

MEASUREMENT	SMALL	MEDIUM	LARGE	X-LARGE
WIDTH <i>side to side</i>	18"	24"	24"	24"
LENGTH /HEIGHT <i>top to bottom</i>	18"	24"	36"	48"

**GONADAL GUARDS**

MEASUREMENT	SMALL	MEDIUM	LARGE
WIDTH <i>side to side</i>	4"	6"	7.5"
LENGTH <i>top to bottom</i>	6"	7"	9"
STRAP <i>edge to edge</i>	8" - 61"	8" - 61"	8" - 61"

**BREAST GUARDS**

MEASUREMENT	Child	Adult
WIDTH <i>side to side</i>	5"	6"
LENGTH <i>top to bottom</i>	6"	7"
STRAP <i>edge to edge</i>	8" - 36"	8" - 36"

**DENTAL DRAPES**

MEASUREMENT	CHILD	ADULT
WIDTH <i>widest point</i>	18"	23.5"
LENGTH <i>from neckline</i>	23.5"	28"
BACK FLAPS <i>from shoulders</i>	7.5"	11.5"

**PANORAMIC DENTAL DRAPES**

MEASUREMENT	ADULT
FRONT WIDTH <i>widest point</i>	24"
FRONT LENGTH <i>from neckline</i>	26"
BACK WIDTH <i>widest point</i>	25"
BACK WIDTH <i>from neckline</i>	15.5"

**SCOLIO GUARDS**

MEASUREMENT	STANDARD
WIDTH <i>side to side</i>	8"
LENGTH <i>top to bottom</i>	26"
NECK STRAP <i>edge to edge</i>	8" - 36"
WAIST STRAP <i>edge to edge</i>	49"

**PREGNANCY FRONT APRONS**

MEASUREMENT	STANDARD	FULL
WIDTH <i>at chest</i>	30"	34"
LENGTH <i>from shoulders</i>	43"	45"

**PREGNANCY HALF APRONS**

MEASUREMENT	STANDARD
WIDTH <i>at hips</i>	27"
LENGTH <i>top to bottom</i>	27"

**LAP GUARD HALF APRONS**

MEASUREMENT	PETITE	SMALL	MEDIUM	LARGE	X-LARGE	2XL
WIDTH <i>at widest point</i>	10"	12"	20.5"	24"	24"	36"
LENGTH <i>top to bottom</i>	10"	12"	16"	18.5"	18.5"	19.5"