

# GONAD/OVARIAN LEAD SHIELD

## PRODUCT INFORMATION


[PHILLIPS-SAFETY.COM](http://PHILLIPS-SAFETY.COM)

COPYRIGHT ©2025. ALL RIGHTS RESERVED



**PHILLIPS SAFETY**  
PRODUCTS INC.

**PRODUCT INFORMATION**  
GONAD/OVARIAN LEAD SHIELD



The Gonad/Ovarian Lead Shield offers a unisex design with easy-to-use buckles and adjustable straps for a secure and comfortable fit. Providing essential protection for the gonadal and ovarian regions during imaging or radiation, these shields are suitable for patients of all genders. Plus, these gonad/ovarian shields have 0.50mm Pb lead equivalency protection level. Available in blue and black, these shields come in a convenient set of three.

**SPECIFICATION**

---

**COLORS:** Black, Blue

**QUANTITY:** Set of 3


**FEATURES:** Adjustable Strap, Buckle


**PROTECTION LEVEL:** 0.50mm Pb lead equivalency

# CONTACT

Should you need any further information,  
please do not hesitate to contact us.

123 Lincoln Blvd, Middlesex, NJ 08846, USA 

+1 866 575 1307 

service@phillips-safety.com 

www.phillips-safety.com 