


# NICU GONAD SHIELDS

PRODUCT INFORMATION



**PHILLIPS SAFETY**  
PRODUCTS INC.

**PRODUCT INFORMATION**  
NICU GONAD SHIELDS



The NICU Gonad Shields protect the reproductive organs of NICU babies during X-ray procedures. These Triangles reduce radiation exposure by 98.7% at 50 kV, and they are both lead-free and latex-free. Designed with infection control as a priority, the Radiation Protection Gonad Shields are disposable and come in a pack of 5. This innovative NICU protective solution significantly lowers the risk of cross-contamination that could occur when using the same lead protection for multiple babies.

**SPECIFICATION**

---

**Core Material: Lead Free**

**Features: Disposable**

**Protection Level: .25mm Pb lead equivalency**

**Size: Small**

**Quantity: Pack of 5**