

Six-Leg, 2-Hook Infusion Pump Stand



MODEL: IV-33

- 2-hook top
- Large diameter chrome plated steel pole
- Oversized base for stability
- 6-leg epoxy coated base
- Single knob height adjustment
- Detachable ram's horn hooks with plastic tips
- 3" rubber wheel, ball bearing casters
- See available options on next page
- 45 lbs. load capacity under normal use


Shipping Weight: 28 lbs.

SPECIFICATIONS

Base Diameter—26"

Height Range: 56"—99"

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

OPTIONS:



IV-46 Pump Tray Support



IV-42 Infusion Pump Frame



IV-50 Pole Outlet Strip




IV-41 Universal Clamp



Heavy Duty Black Epoxy Wire Basket
 IV-51B – 14" x 8" x 8.5"
 IV-52B – 12" x 6" x 6.5"



Heavy Duty Stainless Steel Wire Basket
 IV-51S – 14" x 8" x 8.5"
 IV-52S – 12" x 6" x 6.5"



IV-57 Single Ram's Horn



IV-44 Oxygen Cylinder Holder



IV-56 Drainage Hook



IV-43 Steering Handle



IV-55 Extra Hooks and Organizer Hooks



IV-48 2 lb. Base Stability Weight / V-485 2 lb. Base Stability Weight