

MAMMOGRAPHY/BIOPSY CHAIR



PERFECT CHAIR FOR RADIOLOGY AND MAMMOGRAPHY

MBC000ST MAMMOGRAPHY / BIOPSY

INTENDED USE:

The Hausted® Mammography/Biopsy Chair is ideal for use as a procedural chair or treatment table during mammography/biopsy procedures.



MBC000ST - Product Info

VERSATILITY:

- Its unique layout allows for a full upright position to maneuver between the x-ray machine and the wall
- Self storing side rails allow the chair to maneuver into tight spaces
- Smallest footprint of the Hausted® chairs
- Can be used for per-induction, transport, treatment and recovery
- The pneumatic assisted back rest can adjust for supine to 90° for procedural positioning
- 5 inch Easy Roll casters facilitate easy maneuverability
- Stability Legs added for support during supine position.

PATIENT COMFORT:

- 4-way stretch urethane covered HD foam surface for improved comfort, cleaning and durability
- Leg section and height of the footrest are adjustable, maximizing lower body comfort and stability

| Specifications: | MBC000ST | |
|---|----------|--------|
| | inch | cm† |
| Overall Width: | 27 | 68.58 |
| Overall Length: <i>supine</i> | 75 | 190.5 |
| Overall Length: <i>upright</i> | 42.5 | 107.95 |
| Overall Height: <i>range - minimum</i> | 23.5 | 59.69 |
| Overall Height: <i>range - maximum</i> | 31.5 | 80 |
| <i>patient surface</i> Width: | 22 | 55.88 |
| <i>patient surface</i> Length: | 70.5 | 179.07 |
| Side Rail Length: | 26.74 | 67.92 |
| Side Rail Height: | 10.75 | 27.31 |
| <i>easy rolling</i> Tente® Casters: | 5 | 13 |
| degree of movement ° | | |
| <i>pneumatic assist</i> Backrest: | 0 - 90 | |
| <i>standard</i> Trendelenburg: | 10 | |
| Weight Capacity: EVENLY DISTRIBUTED | lb | kg† |
| | 325 | 147.42 |
| Shipping Weight: | 297 | 134.72 |

† standard rounding method used to show metric values as whole numbers

Product Hyperlinks Enabled



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com



© 2014, GF Health Products, Inc. All Rights Reserved.
Hausted is a registered trademark of GF Health Products, Inc.

graham
field

GF1400052 Rev E22