



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

OptimAL Line, Aluminum Six Drawer MRI-Compatible Emergency Cart, Breakaway Lock, Specialty Package #MR6B-EMG



Specifications

Construction:

- Non-ferrous painted aluminum cabinet and drawers

Features:

- Designed for compatibility in MRI suites
- Tested for safe use in MRI suites by independent testing authority (copies of report available)
- Cabinet painted gloss white with MRI identification sticker applied
- Multiple drawer colors available
- MRI-compatible accessories available
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; full swivel – one with brake, one directional
- Painted aluminum breakaway locking system
- Breakaway lock seals, numbered, yellow, pkg-100
- Corner bumpers, protecting the cart
- Painted aluminum pull-out shelf
- Replaceable plastic top with integrated push handles
- Cart is delivered fully assembled with pre-threaded accessory mounting holes
- Twelve year warranty
- Latex free
- Durable, hard-baked powder coat finish

Assembled Cart Dimensions

- 45.75" H x 34.5" W x 23.75" D
(116.2 cm H x 87.63 cm W x 60.3 cm D)
- Cart net weight: 126 pounds (57.3 kg)
- Approx. shipping weight: 158 pounds (71.8 kg)
- Approx. shipping dims: 57" H x 37" W x 26" D

Internal Drawer Dimensions:

- 23.9" W x 17" D (60.8 cm x 43.2 cm)

Available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Emergency Specialty Package Includes MRI-Compatible Accessories:

- Stainless steel IV pole
- Plastic cardiac board and stainless steel mounting brackets
- Stainless steel oxygen tank holder – D/E size cylinders
- Aluminum utility hooks (4)
- Plastic full drawer trays with adjustable dividers (2)

Drawer Configuration:

- (1) - 4" (10.2 cm)
- (2) - 3.25" (8.3 cm)
- (2) - 6.5" (16.5 cm)
- (1) - 9.75" (24.8 cm)