



Papoose Boards Instructions for Use and Cleaning

INSTRUCTIONS FOR USE:

- Three sets of opposing canvas flaps with hook and loop fasteners fold over the patient's body to provide a snug fit
- Any of the pairs of flaps may be opened while the patient is still held securely by remaining flaps and/or head & arm straps

CLEANING:

- Wipe the board and flap sets with liquid disinfectant or soap and water
- Flap sets and arm straps can also be machine washed

STORAGE:

- When not in use, boards should be stored flat. Normal room temperature and relative humidity is ideal.

MAINTENANCE:

- No end user maintenance required or recommended other than cleaning.

Made in USA

Phillips Safety Products
123 Lincoln Blvd., Middlesex NJ 08846

service@phillips-safety.com
1-732-356-1493