

# Gonad Shields

## PGS-X1



- Set of three or available individually: Infant, Youth and Adult
- 3/4" tie strap
- Easy fastening around the waist
- .50mm LE
- Available in our Classic Nylon Colors

Model Number	Description
PGS-S1	3" x 5 " Infant Size Gonad Shield
PGS-M1	6" x 7 " Youth Size Gonad Shield
PGS-L1	7" x 9 " Adult Size Gonad Shield
PGS-X1	Set of Gonad Shield (Infant, Youth and Adult sizes included)