

# MATERIAL SAFETY DATA SHEET

## 1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

Manufacturer's Name: Jessop Precision Products, Inc. DBA: Suremark  
Address: 4645 Industrial Street, Suite 2C  
City, State, and Zip: Simi Valley, CA 93063  
Emergency Telephone Number: 805-579-9672  
Other Information Calls: 800-406-3399  
Date Prepared: March 2, 2015  
Product Name: Suremark Precision Skin Marking Systems  
To Be Used As: Surface Treatment Markers

## 2. COMPOSITION/INFORMATION ON INGREDIENTS

### 2.1 The following Suremark products:

SL-10, SL-12, SL-15, SL-18, SL-20, SL-25, SL-30, SL-40, SL-50, BL-15, BL-20, SL-15T,  
SL-20T, SL-25T, SL-CR1, SL-W3, SL-W4, BL-W4, SL-M2, SL-M3, SL-W8, SL-W12,  
SL-W15, SL-ARROW, SF-10, SL-LEFT, SL-RIGHT, CC-15, SL-15RT, SL-20RT, SL-25RT

Contains the following toxic chemicals subject to the reporting requirements of Section 313 of the Emergency Planning and Community Right-To-Know Act of 1986 and of 40 CFR 372:

Lead 99.99% by weight CAS# 7439-92-1

### 2.2 All other Suremark products not listed above:

Individual components of this product do not CAS# N/A  
Contribute to any known hazards  
OSHA HAZARDOUS COMPONENTS (20 CFR1910.1200): NONE

## 3. HAZARDS IDENTIFICATION

POTENTIAL HEALTH EFFECTS: NONE

ACUTE: No significant toxic hazard. When burned may produce toxic fumes.

CHRONIC: No significant toxic hazard.

## 4. FIRST AID MEASURES

Emergency and First Aid Procedures:

ROUTES OF ENTRY –

1. Inhalation: N/A
2. Eyes: N/A
3. Skin: N/A

4. Ingestion N/A

## 5. FIRE FIGHTING MEASURES

Flash Point: N/A

Method Used: N/A

Auto-Ignition Temperature: N/A

Flammable Limits in Air % by Volume:

LEL: N/A

UEL: N/A

Extinguishing Media: Use multi-rated ABC extinguisher or water

Special Fire Fighting Procedures: In case of fire, use self-contained breathing apparatus and full protective equipment – all fires liberate toxic gases.

Unusual Fire and Explosion Hazards: None

## 6. ACCIDENTAL RELEASE MEASURES

N/A

## 7. HANDLING AND STORAGE

No special handling required, other than storing away from excessive heat sources.

## 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Respiratory Protection: Material cannot be inhaled under normal conditions of Use. No long-term or irreversible adverse health effects are known.

Skin Protection: Prolonged direct skin contact may cause minor, temporary Irritation and/or dermatitis. No long-term or irreversible Adverse health effects are known.

WARNING! The following products: SL-Left, SL-Right, and TM-25 contain NATURAL RUBBER LATEX in the adhesive. People with LATEX ALLERGIES should avoid contact.

Ingestion: Material is not likely to be ingested under normal conditions of use. Minor throat irritation may occur due to mechanical abrasion of soft tissues if material is swallowed. No long-term or irreversible adverse health effects are known.

Eye Protection: None required

Exposure Guidelines: None established

Engineering Controls: None required

## 9. PHYSICAL AND CHEMICAL PROPERTIES

Boiling Point:	N/A	Specific Gravity (H 0 = 1):	N/A
Vapor Pressure (mm HG):	N/A	Vapor Density (air =1):	N/A
Solubility in Water:	N/A	Reactivity of water:	N/A
Appearance and Odor:	N/A	Melting Point:	N/A

## 10. STABILITY AND REACTIVITY

STABILITY:	N/A
CONDITIONS TO AVOID:	Fire or excessive heat
INCOMPATIBILITY:	Avoid contact with strong oxidizing agents
HAZARDOUS DECOMPOSITION PRODUCTS:	None known
HAZARDOUS POLYMERIZATION:	Will not occur.

## 11. TOXICOLOGICAL INFORMATION

Carcinogen Status ... NTP: NONE    ... IRAC: NONE    ...OSHA: NONE

## 12. ECOLOGICAL INFORMATION

No data available

## 13. DISPOSAL CONSIDERATIONS

WASTE DISPOSAL:    Dispose of in accordance with local, State and Federal regulations.

## 14. TRANSPORT INFORMATION

TRANSPORTATION AND HAZARDOUS MATERIALS DESCRIPTION: Not a hazardous material for DOT shipping.

## 15. REGULATORY INFORMATION

OSHA STANDARD CFR 29, PART 1910.1200: Not listed; non-hazardous

SARA 313 INFORMATION: This product contains no substances at or above the reporting threshold subject to the reporting requirement of section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 and 40 CFR Part 371, based on available data.

CALIFORNIA PROPOSITION 65: The following statement is made in order to comply with the California Safe Drinking Water and Toxic Enforcement Act of 1986. WARNING: This product

contains a chemical known to the State of California to cause cancer and birth defects or other reproductive harm.

## **16. OTHER INFORMATION**

The information in this document is believed to be correct as of the date issued. HOWEVER, NO WARRANTY OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR ANY OTHER WARRANTY IS EXPRESSED OR IS TO BE IMPLIED REGARDING THE ACCURACY OR COMPLETENESS OF THIS INFORMATION, THE RESULTS TO BE OBTAINED FROM THE USE OF THIS INFORMATION OR THE PRODUCT, THE SAFETY OF THIS PRODUCT, OR THE HAZARDS RELATED TO ITS STORAGE, HANDLING, USE, OR DISPOSAL. This information and product are furnished on the condition that the person receiving them shall make his/her own determination as to the suitability of the product for his particular purpose and on the condition that he/she assumes the risk of the use thereof. Jessop Precision Products, Inc assumes no legal responsibility for loss, damage, or expense arising out of, or in any way connected with, the handling, storage, use, or disposal of the product.