

Bariatric, Chrome, Step Stool with Handrail



MODEL: T-6150


- 1" handrail 36" high
- Durable chrome finish
- Welded reinforced "X" brace
- 7/8" diameter handrail 36" high
- 600 lbs. load capacity under normal use

Shipping weight: 15 lbs.

SPECIFICATIONS

PLATFORM: 14.25" L X 11.25" W X 9"H

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1