




UnderShield™  
BY INFAB

# UnderShield™

## Targeted Axillary Protection



UnderShield™ provides targeted radiation protection for the axillary region, addressing the critical gap left by standard aprons. This innovative under-apron solution ensures comprehensive coverage for healthcare professionals working in high-radiation environments.

- Targeted axillary protection regardless of vest size
- Adjustable and flexible design for customized fit
- Exclusive brrr® cooling technology
- Moisture management and rapid drying
- Worn under existing radiation aprons
- Does not limit range of motion

Model / Size	Chest Measurement
US-S	XS/S 33-36" (83.8–91.4 cm)
US-M	M 37-40" (94.0–101.6 cm)
US-L	L 40-43" (101.6–109.2 cm)
US-XL	XL 44-48" (111.8–121.9 cm)
US-2XL	2XL 49-52" (124.5–132.1 cm)
US-3XL	3XL 52-55" (132.1–139.7 cm)

#### Cleaning Instructions:

Remove side panels - do not wash. Vest is machine washable on delicate cycle with mild detergent. Use cold water and hang to dry.

Model Numbers	US-S, US-M, US-L, US-XL, US-2XL, US-3XL
Sizes	S, M, L, XL, 2XL, 3XL
Inner Material	Exclusive brrr® cooling technology
Coverage	Targeted axillary protection
Application	Worn under existing aprons
Comfort	Cooling, moisture-wicking, rapid drying
Protection	Lead Free (Each panel 0.25 mm LE, 0.50 mm LE where overlapped)



REV002

1040 Avenida Acaso  
Camarillo, CA 93012, USA

tel: 805.987.5255  
fax: 805.482.8424

info@infabcorp.com  
www.infabcorp.com

